

LabTech



MultiVap 6

BASIC CONCENTRATOR

USER MANUAL



Thank you for selecting our MultiVap 6.

We are sure that you will be completely satisfied with the performance of this new unit entering your laboratory. We invite you to read carefully this user manual and to keep it close to the instrument for convenient and fast consulting. For any possible clarification or any request for assistance please contact either our local Representative or:

LabTech Srl

Via Fatebenefratelli, 1/5

24010 Sorisole (BG) Italy

Tel +39 035 576614

Website www.labtechsrl.com

E-mail customer.care@labtechsrl.com

INDEX

1. INTRODUCTION.....	4
Features.....	5
2. SAFETY RULES	6
3. INSTALLATION	7
4. SYSTEM SETTING	10
5. TEMPERATURE CURVE	11
6. SOLVENTES TABLE.....	15
7. SERVICE	16

The information contained in this document may be the object of patent applications by LABTECH.

The possession of this document does not confer any license rights in and to such patents.

The following names are LABTECH trademarks throughout the world:

LABTECH
MultiVap 6

All Reproduction Rights Reserved

1. INTRODUCTION

The MultiVap 6 system is a 6-channel parallel nitrogen evaporator developed by LabTech. It is one of the sample preparation instruments for general-purpose and/or standard laboratories.

Compared with a traditional rotary evaporator, the MultiVap 6 system with 6-channel can process 6 parallel samples simultaneously, thereby increasing the laboratory throughput sample processing capability.

The nitrogen flow rate of each channel can be independently adjusted. The nitrogen purge needle in each channel can be adjusted up and down independently.

The system features a smart temperature controller to set and calibrate the temperature via thermocouples, using a robust resistance to heat the aluminum block that optimizes the temperature uniformity in each channel. The temperature difference is less than 5%.

The sample volume is up to 40 mL per channel. The concentration tube is 40 mL of glass vial with graduation and can be used for different instruments. Without transferring, samples are directly concentrated, thereby reducing transfer loss and improving recovery.

The system is compact (only 300 mm by 170 mm and 120 mm) and light with a weight of only 5 kg and can be easily placed in a fume hood.

Features

Channel number: 6 parallel.

Heating mode: aluminum block dry heating.

Nitrogen purge: purge needle can be adjusted individually up and down as well as the nitrogen flowrate that can be adjusted individually.

Concentration tube: 40 mL graduated glass flat-bottom, 40 mL tail-tube with scale (1mL end-point).

Anti-corrosion material: PTFE or Teflon are used for the solvent exposed parts, can endure the dichloride methonal solvent.

Instrument dimensions: 300 (W) × 170 (H) × 120 (D) mm.

Power: AC 230V, 50/60Hz.

Gas supply: nitrogen 0.1—0.2MPa.

Weight: 5 kg.

Max Temperature: 120°C.

Supplied with each MultiVap 6:

Collection bottle 40 mL: 6.

Power cord 220V 50/60Hz: 1.

Two-way connector O.D4mm/O. D6 mm: 1.

PU tubing O.D4 mm, 2 m: 1.

PU tubing O.D6 mm, 2 m: 1.

User manual: 1.

2. SAFETY RULES

General Information

Please read carefully this user manual before starting to use the instrument and follow its prescriptions with the utmost care. This user manual is part of the delivery, hence must be always kept together with the instrument on its working place.

It is imperative that every person operating with this system has read and fully understood this manual. The non-observance of the instructions contained herein or improper use may involve damages/injuries that are not covered by product liability.

Electrical safety



The instrument has to be used within the rated voltage. Prior to use, please check if the wire is aged. In case of aged wires, please contact the after-sales service for inspection. It is forbidden to disassemble the instrument and to connect internal circuit parts, in order to avoid a short circuit or open circuit.

Fire safety

Numerous reagents are flammable and explosive. When the solvent vapor concentration reaches a certain level, it would be flammable and could cause fire. The instrument should be kept away from the sources of ignition and high temperature places. If there is solvent pungent smell, carefully check whether there is gas or liquid leakage, and turn off the power.



Chemical safety


The instrument is an instrument for organic chemical sample pretreatment. The involved chemical solvents have harmful effects on the human health. Despite the instrument is fully closed and features full vent design, it is recommended to pay attention to the personal safety during the use. Regular check of liquid waste barrels as well as working conditions of the vent fan are required to avoid the risk of leakage caused by corrosion and to avoid the formation of organic solvent vapors affecting operators' health. If there is a fault, please contact the Labtech Service Team.



3. INSTALLATION

Setting up your LabTech MultiVap 6 Just follow the following instructions. Once completed the setup, take some time to explore the features of your system. The user manual provides tips and instructions to help you learn the basics of your MultiVap 6

TOP VIEW




The instrument consists of six channels, the gas flow rate and stainless-steel purge needle height can be independently adjusted in every channel.

The red painted part of the instrument is the heating aluminum block, which can uniformly heat six channels. The temperature is controlled by smart temperature controller, easy to use, easy and accurate to be set and calibrated. The temperature difference is less than 5%.

Concentration tube: 40 mL flat-bottom or tail tube with graduation (optional).

TOP VIEW



1 - Nitrogen tubing: FEP transparent tubing as the gas feeding pipe. FEP is a kind of polytetrafluoroethylene (PTFE), which is good anti-corrosive material for various solvents.

2 - Flanged/screwed joints: PP + fiberglass flange joints with good strength and solvent resistance, for connecting the purge needle and gas tubing.

3 - Purge needle holder: to control the purge needle adjusted up-down and left-right regulation.

4 - Concentration tubes: 40 mL graduated glass flat-bottomed or with tailpipe concentration vial with graduation.

5 - Purge needle: stainless steel purge needle can be adjusted up and down, left and right, to find the best concentrating position.

6 - Main body: MultiVap 6 nitrogen blowing concentrator main sheet metal shell, with insulated pad between heating aluminum block.

7 - Heating aluminum blocks: a special aluminum block used to uniformly heat concentration tubes.

8 - Smart temperature controller: temperature setting and calibration.

9 - Gas flow indicator: unscrew to increase the pressure, screw to lower the pressure.

10 - Gas flow adjustment knob: one knob per channel.

TOP VIEW

- 1 - power switch.
- 2 - three-phase power socket.
- 3 - 4 mm diameter gas inlet.

4. SYSTEM SETTING











- 1- Top display window (showing measuring value PV, parameter name).
- 2 - Bottom display window (showing the setting value SV, alarm code, parameter data, etc.).
- 3 - Setting button (to enter parameter setting, and confirming, etc.).
- 4 - Data moving button (may be used to integrate data).
- 5 - Data reducing button.
- 6 - Data increasing button.
- 7 - LED indicator (OP1/AU1/AU2/RUN showing data output, alarm setting 1,2 and running indicator).



Basic display status:

when the power is on, the display window displays the measurement value PV (top) and setting value SV (bottom). When the measurement value is out of range of the input value (such as thermocouple broken) "orA" is alternately displayed on the top window with measuring upper or lower limit and the instrument will be automatically stopped.

Set the given values:

in the basic display status, if the parameter lock is unlocked, the temperature control value can be modified on the bottom window by pressing ,  or , pressing  to decrease the data and pressing  to increase the data. The decimal point would flash for modifying data. Holding the buttons can quickly increase or decrease the value and speeding up with the decimal point shifted to the right. Pressing  can move the location of the data (cursor) directly. Pressing  or  can modify the value of the flashing position and the operation is very quick.

5. TEMPERATURE CURVE








The temperature is increased quickly in the first 10 minutes, up to 4°C / more or less.
 The minimum RSD between two channels: 0.11%, indicating that the temperature difference between channels is low.
 There is about 5°C difference between the setting value and the real one.




8 SOLVENTS TEST

Parallel test among CH1-CH6 with conducting time and volume by using 5, 10, 20, 30 and 40 mL solvent as test solution, using a graduated cylinder to measure the volume and a stopwatch to measure the time.

Solvents	Boiling point°C	Solvents	Boiling point°C
Methanol	64.8	Cyclohexane	80.7
Acetonitrile	81.1	Methylene chloride	39.8
Ethyl acetate	77.2	Acetone	56.5
Hexane	68.74	Petroleum ether	30-60







6. SOLVENTS TABLE

Solvent	Formula	Molar mass	Evaporation energy	Boiling point	Density in	Vacuum in mbar for
		in g/mol	in J/g	at 1013 mbar	g/cm ³	boiling point at 40 °C
Acetone	CH ₃ C(=O)CH ₃	58.1	553	56	0.790	556
n-amylalcohol, n-pentanol	C ₅ H ₁₂ O	88.1	595	37	0.814	11
Benzene	C ₆ H ₆	78.1	548	80	0.877	236
n-butanol	C ₄ H ₁₀ O	74.1	620	118	0.810	25
tert. butanol (2-methyl-2-propanol)	C ₄ H ₁₀ O	74.1	590	82	0.789	130
Chlorobenzene	C ₆ H ₅ Cl	112.6	377	132	1.106	36
Chloroform	CHCl ₃	119.4	264	62	1.483	474
Cyclohexane	C ₆ H ₁₂	84.0	389	81	0.779	235
Diethylether	C ₄ H ₁₀ O	74.0	389	35	0.714	850
1,2-dichloroethane	C ₂ H ₄ Cl ₂	99.0	335	84	1.235	210
1,2-dichloroethylene (cis)	C ₂ H ₂ Cl ₂	97.0	322	60	1.284	479
1,2-dichloroethylene (trans)	C ₂ H ₂ Cl ₂	97.0	314	48	1.257	751
Diisopropyl ether	C ₆ H ₁₄ O	102.0	318	68	0.724	375
Dioxane	C ₄ H ₈ O ₂	88.1	406	101	1.034	107
DMF (dimethyl-formamide)	C ₂ H ₇ NO	73.1		153	0.949	11
Acetic acid	C ₂ H ₄ O ₂	60.0	695	118	1.049	44
Ethanol	C ₂ H ₆ O	46.0	879	79	0.789	175
Ethylacetate	C ₄ H ₈ O ₂	88.1	394	77	0.900	240
Heptane	C ₇ H ₁₆	100.2	373	98	0.684	120
Hexane	C ₆ H ₁₄	86.2	368	69	0.660	360
Isopropylalcohol	C ₃ H ₈ O	60.1	699	82	0.786	137
Isoamylalcohol (3-methyl-1-butanol)	C ₅ H ₁₂ O	88.1	595	129	0.809	14
Methylethylketone	C ₄ H ₈ O	72.1	473	80	0.805	243
Methanol	CH ₃ O	32.0	1227	65	0.791	337
Methylene chloride, dichloromethane	CH ₂ Cl ₂	84.9	373	40	1.327	850
Pentane	C ₅ H ₁₂	72.1	381	36	0.626	850
n-propylalcohol	C ₃ H ₈ O	60.1	787	97	0.804	67
Pentachloroethane	C ₂ HCl ₅	202.3	201	162	1.680	13
1,1,1,2-tetra-chloroethane	C ₂ HCl ₄	167.9	247	146	1.595	20
Tetrachlorocarbon	CCl ₄	153.8	226	77	1.594	271
1,1,1-trichloroethane	C ₂ HCl ₃	133.4	251	74	1.339	300
Tetra-chloro-ethylene	C ₂ Cl ₄	165.8	234	121	1.623	53
THF (tetrahydrofurane)	C ₄ H ₈ O	72.1		67	0.889	374
Toluene	C ₇ H ₈	92.2	427	111	0.867	77
Trichloroethylene	C ₂ HCl ₃	131.3	264	87	1.464	183
Water	H ₂ O	18.0	2261	100	1.000	72
Xylene (mixture)	C ₈ H ₁₀	106.2	389			25
o-xylene	C ₈ H ₁₀	106.2		144	0.880	
m-xylene	C ₈ H ₁₀	106.2		139	0.864	
p-xylene	C ₈ H ₁₀	106.2		138	0.861	

7. SERVICE

The LABTECH worldwide technical support network consists of highly trained Field Service Engineers, Technical Support Specialists and Service Coordinators who are ready to quickly assist customers with answers and solutions to service needs and application questions.

For any possible clarification or any request for assistance please contact either our local Representative or:

LabTech Srl

Via Fatebenefratelli, 1/5

24010 Sorisole (BG) Italy

Tel +39 035 576614

Website www.labtechsrl.com

E-mail customer.care@labtechsrl.com