

HEATING SOLUTIONS

Hotplates

LabTech *EH* and *EG* hotplates are designed to meet laboratory daily heating requirements. Their unique features include different digital temperature controller and anticorrosive or aluminium heating surface.



EH LINE

- Large heating surface. Fast and uniform heating.
- Microprocessor chip for precise temperature control.
- Constant heating and no electromagnetic interference.
- Special high resistant coatings prevent corrosion and guarantee long life.
- Compact design combined with an elegant look.

- Large heating surface. Fast and uniform heating.
- Microprocessor chip for precise temperature control.
- Constant heating and no electromagnetic interference.
- Special high resistant coatings prevent corrosion and guarantee long life.
- Innovative design with stand-alone temperature control box.



EG LINE

EH LINE TECHNICAL PARAMETERS

Model	EH-20A Plus	EH-20B	EH-20D	EH-35A Plus	EH-35B
Heating Power	2000W				
Power Supply	230V-50/60Hz				
Temp. Control Mode	Switch (digital display)	PID (digital display)	PID (digital display)	Switch (digital display)	PID (digital display)
Temp. Stability	±5 °C	±1 °C	±0.2 °C	±5 °C	±2 °C
Heating Mode	Electrothermal Component				
Temp. Range	Ambient ~ 200°C			Ambient ~ 350 °C	
Material	Aluminium alloy with anticorrosive coating			Aluminium alloy and stainless steel	
Corrosion-resistant	Excellent			Good	
Continuing Working Time	>48 hours				
Heating Area	400×300mm				
Dimensions(L×W×H)	400×300×120mm			400×300×150mm	
Weight	18Kg				

EG LINE TECHNICAL PARAMETERS

Model	EG-20A Plus	EG-20B	EG-35A Plus	EG-35B
Heating Power	3000W			
Power Supply	230V-50/60Hz			
Temp. Control Mode	Switch (digital display)	PID (digital display)	Switch (digital display)	PID (digital display)
Temp. Stability	±5 °C	±2 °C	±5 °C	±2 °C
Heating Mode	Electrothermal Component			
Temp. Range	Ambient ~ 200°C		Ambient ~ 350°C	
Material	Aluminium alloy with anticorrosive coating		Aluminium alloy and stainless steel	
Corrosion-resistant	Excellent		Good	
Continuing Working Time	>48 hours			
Heating Area	600×400mm			
Dimension(L×W×H)	600×400×110mm			
Weight	30Kg			

For any question on **LabTech** Hotplates, please write to info@labtechsrl.com

LabTech Srl

Via Fatebenefratelli 1/5
24010 Sorisole (BG), Italy
Tel. +39 035 576614
Fax +39 035 4729414

Tel.: +39 035 576614
Fax: +39 035 4729414
E-mail: info@labtechsrl.com
Web site: www.labtechsrl.com