






## Accessory List

Rotation Flask				
	P/N	Volume	Joint Specification	
	50103-0050	50mL	29/32	
	50103-0100	100mL	29/32	
	50103-0250	250mL	29/32	
	50103-0500	500 mL	29/32	
	50103-1000	1000 mL	29/32	
	50103-2000	2000 mL	29/32	
	50103-3000	3000 mL	29/32	
	50102-0100	100 mL	24/29	
50102-0250	250 mL	24/29		
Heart-shaped Flask				
	P/N	Volume	Joint Specification	
	50102-0050J	50 mL	24/29	
	50102-0100J	100 mL	24/29	
50101-0050J	50 mL	19/16		
Conversion Adapter				
	P/N	Joint Specification		
	50103-29-24	29/32	24/29	
	50103-29-19	29/32	19/26	
50103-24-19	24/29	19/26		
Collection Flask				
	P/N	Volume	Joint Specification	
	50103-C0	500mL	35/20	
	50103-C1	1000mL	35/20	
50103-C2	2000mL	35/20		
Anti-splashing Connector				
	P/N	Volume	Joint Specification	
	50103-0001	100mL	29/32	24/29