


ROTARY EVAPORATOR ACCESSORIES

1. Vertical Condenser / Diagonal Condenser / Dry Ice Condenser



ROTARY EVAPORATORS ACCESSORIES

| Item | Part number | Description |
|------|-------------|--------------------|
| 1 | 50103-V | Vertical Condenser |
| 1 | 50103-L | Diagonal Condenser |
| 1 | 50103-I | Dry Ice Condenser |
| 2 | 50103-NC | Hose Nipple |
| 3 | 50103-A | Stop Cock |
| 4 | 50103-CS | Ball Joint Clamp |
| 5 | 50103-SR3 | Lock Handwheel |
| 6 | 50103-SR9 | Lock Spring |

| Item | Part number | Description |
|------|-------------|----------------------------|
| 7 | 50103-SR1 | Seal Ring Set |
| 8 | 50103-NS | Lock and Push-off nuts Set |
| 9 | 50103-GA | Vapor Tube |
| 10 | 50103-C1 | Clip 19 mm |
| 10 | 50103-C2 | Clip 24 mm |
| 10 | 50103-C3 | Clip 29 mm |
| 11 | 50103-1000 | Evaporating Flask 1L |
| 12 | 50103-C1 | Receiving Flask 1L |

Plastic coating for glassware is available upon request

Evaporation Flask

| P/N | Volume | Joint Size |
|------------|---------|------------|
| 50103-0050 | 50 mL | 29/32 |
| 50103-0100 | 100 mL | 29/32 |
| 50103-0250 | 250 mL | 29/32 |
| 50103-0500 | 500 mL | 29/32 |
| 50103-1000 | 1000 mL | 29/32 |
| 50103-2000 | 2000 mL | 29/32 |
| 50103-3000 | 3000 mL | 29/32 |
| 50103-0100 | 100 mL | 24/29 |
| 50103-0250 | 250 mL | 24/29 |



Trap Ball

| P/N | Volume | Joint Size |
|------------|--------|----------------|
| 50103-0001 | 100 mL | 29/32 24/29 |



Heart-shaped Flask

| P/N | Volume | Joint Size |
|-------------|--------|------------|
| 50102-0050J | 50 mL | 24/29 |
| 50102-0100J | 100 mL | 24/29 |
| 50101-0050J | 50 mL | 19/16 |



Receiving Flask

| P/N | Volume | Joint Size |
|----------|---------|------------|
| 50103-C0 | 500 mL | 35/20 |
| 50103-C1 | 1000 mL | 35/20 |
| 50103-C2 | 2000 mL | 35/20 |



Adaptor

| P/N | Joint Size | Joint Size |
|-------------|------------|------------|
| 50103-29-24 | 29/32 | 24/29 |
| 50103-29-19 | 29/32 | 19/26 |
| 50103-24-19 | 24/29 | 19/26 |



Spider with 5 flasks

| P/N | Volume | Joint Size |
|-------------|----------------|----------------|
| 50104-29-24 | 50 mL 5 pcs | 29/32 24/29 |

