

COVID 19 FIGHT SERUM TEST

CONFIGURATIONS BASED ON 5/7 mL ø 13 mm CAPPED TUBES

SMALL CAPACITY: 20 PLACES

1 x CENTRIFUGE 416 RPM MAX 4.000, RCF MAX 2.826 OR 624R RPM MAX 6.000, RCF MAX 3.960
1 x SWING OUT ROTOR GRS-S-100-4 4 x 100 mL RPM MAX 4.000, RCF MAX 2.774 ON 416
1 x SWING OUT ROTOR GRS-S-100-4 4 x 100 mL RPM MAX 4.200, RCF MAX 3.059 ON 624R
4 x BUCKETS GLB-b100-100 1 x 100 mL
4 x ADAPTERS GAM-7-5(b100) EACH 5 PLACES FOR A TOTAL OF 20
Rotor for MTP and DWP available.

MEDIUM CAPACITY: 32 PLACES

1 x CENTRIFUGE 1248 RPM MAX 12.000, RCF MAX 15.520 OR 1248R RPM MAX 12.000, RCF MAX 16.582
1 x SWING OUT ROTOR GLB-G-250-4 4 x 250 mL RPM MAX 5.000, RCF MAX 5.394
4 x BUCKETS GLB-250-250 1 x 250 mL
4 x ADAPTERS GAM-7-8(250) EACH 8 PLACES FOR A TOTAL OF 32
Rotor for MTP and DWP available. Fixed angle, high speed rotors available.

MEDIUM / LARGE CAPACITY: 48 PLACES

1 x CENTRIFUGE 1248 RPM MAX 12.000, RCF MAX 15.520 OR 1248R RPM MAX 12.000, RCF MAX 16.582
1 x ROTOR GRS-G-r250-4 4 x 250 mL RPM MAX 4.000, RCF MAX 3.134
4 x BUCKETS GLB-r250-r250 1 x 250 mL
4 x ADAPTERS GAM-5-12(r250) EACH 12 PLACES FOR A TOTAL OF 48
Holder for MTP and DWP available. Fixed angle, high speed rotors available.

LARGE CAPACITY: 120 PLACES

1 x CENTRIFUGE 1580 RPM MAX 15.000, RCF MAX 24.249 OR 1580R RPM MAX 15.000, RCF MAX 25.910
1 x ROTOR GRS-G-r750-4 4 x 750 mL RPM MAX 3.600, RCF MAX 2.898
4 x BUCKETS GLB-r750-r750 1 x 750 mL
4 x ADAPTERS GAM-5-30(r750) EACH 30 PLACES FOR A TOTAL OF 120
Holder for MTP and DWP available. Fixed angle, high speed rotors available.

IMPRESSIVE CAPACITY: 192 PLACES !!

1 x CENTRIFUGE 1696R RPM MAX 16.000, RCF MAX 27.590
1 x ROTOR GRS-B-1000-4 4 x 1.000 mL RPM MAX 4.000, RCF MAX 2.939
4 x BUCKETS GLB-h1000-1000 1 x 1000 mL
4 x ADAPTERS GAM-5-48(h1000) EACH 48 PLACES FOR A TOTAL OF 192
Holder for MTP and DWP available. Fixed angle, high speed rotors available.