

WATER CHILLERS COOLING SOLUTIONS

LabTech

LARGE COOLING CAPACITY WATER CHILLER



COOLING SOLUTIONS BY LABTECH

Even if water is one of the most precious and scarce resources it is still largely used as coolant in many laboratories.



LabTech offers a complete line of recirculating water chillers to match all laboratory needs. The LabTech's state-of-the-art water chillers ensure accurate and constant cooling conditions of instruments even in the harsh laboratory conditions.

LabTech follows "Green Lab Conditions" regulations by using eco-friendly materials by adapting new technologies and innovations to maintain top reliability level worldwide. All production steps are focused to offer high quality and customized solutions to meet any requirement.

The LabTech water chiller line is specially designed for analytical, medical and industrial use providing accurate control of the temperature.

TECHNICAL SPECIFICATIONS

H150-10000N	
Temperature stability (°C)	+8 ~ +35
Temperature stability (°C)	± 0.1
Temperature control mode	Microprocessor
Cooling Capacity (W)	10000
Pump flow (L/min)	28 @40psi
Pump power (W)	600
Pump type	Centrifugal
Heat exchange model	Brazed Plate
Dimension (LxWxH cm)	71x74x117
Weight (Kg)	145
Electric power supply	400V 50Hz
Order information	LWGCA9

MAIN APPLICATIONS

- X-ray
- Lasers
- Industrial
- Liquid Nitrogen Generators
- Electron microscope
- Centralized Cooling Circuits

LabTech

WWW.LABTECHSRL.COM

LabTech, your reliable laboratory solution provider.

LabTech is a leading manufacturing company with skill and passion devoted to provide advanced laboratory solutions to the analytical community. We are located in Sorisole, northern Italy, with facilities in USA and China. Incorporated with global branches, LabTech is organized into R&D, production, marketing & sales and customer support. Major LabTech products include organic/inorganic sample preparation systems, laboratory cooling/heating/temperature control equipments, analytical instruments and gas generators. With knowledgeable scientists and experienced engineers, our R&D team has developed and released a series of new lab equipments and instruments dedicated to the modern analytical laboratories. With all of these accomplishments and our endless research and development efforts, we are proud to be your professional partner.



LabTech Srl

Via Fatebenefratelli, 1/5 - 24010 Sorisole (BG) - Italy
Tel: +39 035 576614 - Fax: +39 035 4729414
www.labtechsrl.com - email: marketing@labtechsrl.com